



PORTABLE AUTOCLAVE

We are engaged in offering high quality Autoclave Sterilizers, which are manufacture using stainless steel. Steam autoclave is a device that is used for heating the substances above their boiling point and also used to sterilize the surgical instruments. These are extensively used for the effective cleaning of surgical and laboratory instruments for their safety and hygienic from pesticides.

SALENT FEATURES



- Available in circular design (Optional : Rectangular design)
- Versatile usage
- Ergonomic & Sturdy Design
- Energy Efficient
- Low Maintenance & long life
- Autoclave S. S. body
- Vertical wing nut type seamless
- Fitted with pressure gauge
- 2 Spring Loaded Safety Valve
- Vacuum breaker for safety
- Hydraulically tested up to 40 psi

SPECIFICATION

Model	PA-100	PA-110	PA-120
Sterilizing Temp	121 °C	121 °C	121 °C
Sterilizing Time	15-20 minutes	15-20 minutes	15-20 minutes
Heating Element	1500 W, 230 V AC	2000 W, 230 V AC	2000 W, 230 V AC
Capacity (Liters)	6	12	17
Chamber dimension	200 x 200 mm	250 x 250 mm	280 x 280 mm
Working pressure	1.1 to 1.2 Kg/cm2 (15-18PSI)	1.1 to 1.2 Kg/cm2 (15-18PSI)	1.1 to 1.2 Kg/cm2 (15-18PSI)

Optional:

- Digital PID temp indicator cum controller along with timer
- Pressure Control Switch can be fitted.
- Digital indicator cum controller can be fitted.
- Digital timer can be fitted.
- Water cut off devices can be fitted.

Note: Customized Solution Available On Request