



REAGENT REFRIGERATOR



SAILENT FEATURES

- Inner chamber made of stainless steel-304 and outer chamber made of mild steel duly powder coated.
 - PUF thick insulation to ensure temperature stability & reduced energy consumption
 - Keyed door lock for improved security with magnetic door gasket.
 - SS/plastic coated adjustable trays are provided for storage inside (varies on size of refrigerator and demand). Baskets can also be provided as per requirement
 - Micro - processor controls the internal temperature between 2°C to 8°C with an accuracy of $\pm 1^\circ\text{C}$
 - Specially designed controller for audio visual deviation alarm from the preset temperature, door open alarm and power failure alarm and 4 x 7 segment LED display.
 - Hot line around the mouth of cabinet to prevent moisture condensation in humid weather
 - The cooling is effected by a hermetically sealed CFC FREE ECO FRIENDLY COMPRESSOR using refrigerant R134a.
 - 8 Watt Fluorescent lamp / CFL illumination
- It works on 230V AC.
 - Designed air flow system ensuring uniform temperature distribution
 - Unit mounted on castor wheels for easy of movement
 - Frost free design

SPECIFICATION

Model	SAII-RF-60	SAII-RF-110	SAII-RF-180	SAII-RF-250	SAII-RF-340	SAII-RF-500
Capacity (Liters)	60	110	180	250	340	500
Internal dimension (WxDxH cms)	30 x 30 x 65	40 x 40 x 70	50 x 45 x 80	55 x 51 x 90	60 x 57 x 100	65 x 65 x 120
External dimension (WxDxH cms)	50 x 50 x 115	60 x 60 x 120	70 x 65 x 135	75 x 71 x 160	80 x 77 x 170	85 x 85 x 190
Temperature range	2 to 8°C					

Optional:

- Circular chart temperature recorder
- Software to interface freezer directly to computer for storing data & print
- In built voltage stabilizer