



SEROLOGICAL WATER BATH (DOUBLE WALL)

Salient Features:

- Double walled water bath
- Inside temperature from to Ambient + 5°C to 95°C can be set freely
- Outer body is powder coated CRC steel sheet
- Inner chamber made of SS 304.
- High grade mineral wool insulation to avoid thermal loss
- Thermostat controlled with an accuracy of $\pm 1^{\circ}\text{C}$
- Neon lamp indicator main switch
- Top opening.
- Cover is provided at the top with perforated shelf



Specifications:

Model	Size	Capacity	Temperature range	Weight (Kg)	Heating load
SAII-SWB-100	12" x 6" x 6"	6 Liters	Ambient +5°C to +95°C	7	1 kW
SAII SWB-110	12" x 10" x 7"	13 Liters	Ambient +5°C to +95°C	10	1.5 kW
SAII SWB-120	14" x 10" x 7"	16 Liters	Ambient +5°C to +95°C	12	1.5 kW
SAII SWB-130	18" x 10" x 8"	23 Liters	Ambient +5°C to +95°C	14	2 kW
SAII SWB-140	24" x 10" x 8"	32 Liters	Ambient +5°C to +95°C	12	2 kW
SAII SWB-150	24" x 14" x 12"	66 Liters	Ambient +5°C to +95°C	20	3 kW

Optional:

- SS basket
- Stirrer can be fitted at an extra cost
- SS outer body
- PID temperature indicator cum controller
- Automatic digital timer
- RS 232 port and communication software to interface with PC

NOTE: CUSTOMIZED SOLUTION ALSO AVAILABLE ON REQUEST